

# HidraWear

REDEFINING PATIENT LED WOUND CARE

**HidraWear** is the worlds first **HS specific** wound care product designed to make managing HS drainage quick and easy, and completely adhesive free. It is a unique Wearable Wound Care system for the management of wounds in difficult to dress areas, such as the under-arm region, or buttocks, groin & upper thigh.

Our **NEW** dressing retention system combines reusable garments with a super absorbent dressing requiring no adhesive and is secured using our patented HidraHex technology.

Each garment utilizes our patented HidraHex technology to enable patients to easily change their dressing day or night

## HidraWear Garments Incorporate:



### HidraHex

The hexagon pattern is designed to minimize leaks and strike through while remaining fully breathable & promoting wicking.



### SeamSense

All HidraWear garments incorporate outward facing seams in traditional areas of high friction. This ensures a super soft skin contacting surface where it matters most, reducing risk of irritation to provide smooth movement.



### Easy On

HidraWear Garments are specifically designed to be easy to get on and off, with the women's underarm crop top incorporating a front closing clasp



### AeroSilk

All HidraWear garments are made of a premium micro modal and elastane blend fabric that is soft and breathable.



### TrueFit


Adjustable straps in the women's underarm crop top enable a tailored fit.



**Retention Dressing System includes a washable garment and non adhesive super absorbent dressings**

### HidraWear Dressing

3 x 5 inch super absorbent dressing that locks away moisture from the wound while protecting clothing and bed-linen from exudate strike through. **HidraHex Hook & Loop** retention technology secures the dressing in place without use of any adhesives.



HidraWear can now be prescribed to your patients, with their monthly supply delivered directly to their home

Simply fill in the order from on page 3 & email or fax along with your patient info page (See right for filled in example)

## How It Works:



Delivered to patient

<b>Clinic Name</b> ABC Dermatology <small>Clinic Phone 123-456-7890 Clinic Fax 123-456-7890</small>																																																																																												
<b>Patient Name</b> Jane Doe		<b>FAX # (866) 220-8963</b> <small>Fax or email* with Patient Information Sheet orders@rotech.com Phone # (888) 711-2014</small>																																																																																										
<b>Order Date:</b> 11-2-2022		<small>AE Name</small> _____ <small>AE ID #</small> _____																																																																																										
<small>Have patient's wound/s ever been debrided?</small> <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO		<small>Is patient currently seen by Home Health?</small> <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO																																																																																										
<table border="1"> <tr> <td><small>ICD-10</small></td> <td><small>L73.2-Hidradentis suppurativa</small></td> <td><small>L98.499-Non-pressure chronic ulcer of skin</small></td> <td><small>L98.419-Non-pressure chronic ulcer of skin / Buttocks</small></td> </tr> </table>				<small>ICD-10</small>	<small>L73.2-Hidradentis suppurativa</small>	<small>L98.499-Non-pressure chronic ulcer of skin</small>	<small>L98.419-Non-pressure chronic ulcer of skin / Buttocks</small>																																																																																					
<small>ICD-10</small>	<small>L73.2-Hidradentis suppurativa</small>	<small>L98.499-Non-pressure chronic ulcer of skin</small>	<small>L98.419-Non-pressure chronic ulcer of skin / Buttocks</small>																																																																																									
<table border="1"> <thead> <tr> <th>ICD-10 (See ICD-10 legend above)</th> <th>Wound Thickness (HidraWear dressing requires Full)</th> <th>LOCATION</th> <th>Wound Size (cm) (L x W x D)</th> <th>Drainage (Exudate) HidraWear Dressing requires M - H</th> </tr> </thead> <tbody> <tr> <td>Wound 1: L73.2 + 98.499</td> <td><input type="checkbox"/> Partial <input checked="" type="checkbox"/> Full</td> <td>Right Armpit</td> <td>2.5 x 1.2 x 0.2</td> <td>N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input checked="" type="checkbox"/></td> </tr> <tr> <td>Wound 2: L73.2 + 98.419</td> <td><input type="checkbox"/> Partial <input checked="" type="checkbox"/> Full</td> <td>Left Buttocks</td> <td>4.6 x 2.8 x 0.2</td> <td>N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input checked="" type="checkbox"/></td> </tr> <tr> <td>Wound 3:</td> <td><input type="checkbox"/> Partial <input type="checkbox"/> Full</td> <td></td> <td></td> <td>N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/></td> </tr> <tr> <td>Wound 4:</td> <td><input type="checkbox"/> Partial <input type="checkbox"/> Full</td> <td></td> <td></td> <td>N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/></td> </tr> <tr> <td>Wound 5:</td> <td><input type="checkbox"/> Partial <input type="checkbox"/> Full</td> <td></td> <td></td> <td>N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/></td> </tr> <tr> <td>Wound 6:</td> <td><input type="checkbox"/> Partial <input type="checkbox"/> Full</td> <td></td> <td></td> <td>N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/></td> </tr> <tr> <td>Wound 7:</td> <td><input type="checkbox"/> Partial <input type="checkbox"/> Full</td> <td></td> <td></td> <td>N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/></td> </tr> <tr> <td>Wound 8:</td> <td><input type="checkbox"/> Partial <input type="checkbox"/> Full</td> <td></td> <td></td> <td>N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/></td> </tr> </tbody> </table>				ICD-10 (See ICD-10 legend above)	Wound Thickness (HidraWear dressing requires Full)	LOCATION	Wound Size (cm) (L x W x D)	Drainage (Exudate) HidraWear Dressing requires M - H	Wound 1: L73.2 + 98.499	<input type="checkbox"/> Partial <input checked="" type="checkbox"/> Full	Right Armpit	2.5 x 1.2 x 0.2	N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input checked="" type="checkbox"/>	Wound 2: L73.2 + 98.419	<input type="checkbox"/> Partial <input checked="" type="checkbox"/> Full	Left Buttocks	4.6 x 2.8 x 0.2	N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input checked="" type="checkbox"/>	Wound 3:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>	Wound 4:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>	Wound 5:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>	Wound 6:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>	Wound 7:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>	Wound 8:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>																																												
ICD-10 (See ICD-10 legend above)	Wound Thickness (HidraWear dressing requires Full)	LOCATION	Wound Size (cm) (L x W x D)	Drainage (Exudate) HidraWear Dressing requires M - H																																																																																								
Wound 1: L73.2 + 98.499	<input type="checkbox"/> Partial <input checked="" type="checkbox"/> Full	Right Armpit	2.5 x 1.2 x 0.2	N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input checked="" type="checkbox"/>																																																																																								
Wound 2: L73.2 + 98.419	<input type="checkbox"/> Partial <input checked="" type="checkbox"/> Full	Left Buttocks	4.6 x 2.8 x 0.2	N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input checked="" type="checkbox"/>																																																																																								
Wound 3:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>																																																																																								
Wound 4:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>																																																																																								
Wound 5:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>																																																																																								
Wound 6:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>																																																																																								
Wound 7:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>																																																																																								
Wound 8:	<input type="checkbox"/> Partial <input type="checkbox"/> Full			N <input type="checkbox"/> L <input type="checkbox"/> M <input type="checkbox"/> H <input type="checkbox"/>																																																																																								
<b>Garment Sizing</b> <table border="1"> <tr> <td><b>Crop Top</b></td> <td colspan="5">Bias decision to bust size - If between sizes, downsizing is preferred</td> </tr> <tr> <td></td> <td>Small</td> <td>Medium</td> <td>Large</td> <td>X-Large</td> <td>XX-Large</td> </tr> <tr> <td>Bust (Inch)</td> <td>32-36</td> <td>36-40</td> <td>36-44</td> <td>44-48</td> <td>48-52</td> </tr> <tr> <td>Waist (Inch)</td> <td>24-28</td> <td>28-32</td> <td>32-36</td> <td>36-40</td> <td>40-44</td> </tr> </table>		<b>Crop Top</b>	Bias decision to bust size - If between sizes, downsizing is preferred						Small	Medium	Large	X-Large	XX-Large	Bust (Inch)	32-36	36-40	36-44	44-48	48-52	Waist (Inch)	24-28	28-32	32-36	36-40	40-44	<p>** One Garment provided with first order</p> <p>HidraWear is the first and only Hidradentis Suppurativa (HS) specific wound dressing system, intended for use by people with wounds in difficult to dress areas of the body such as the armpit, buttocks &amp; groin. It is comprised of a reusable washable garment along with ergonomically designed single-use wound pads that come with a patented dressing fastener. HidraWear is clinically proven to assure secure dressing retention without the use of adhesives</p>																																																																		
<b>Crop Top</b>	Bias decision to bust size - If between sizes, downsizing is preferred																																																																																											
	Small	Medium	Large	X-Large	XX-Large																																																																																							
Bust (Inch)	32-36	36-40	36-44	44-48	48-52																																																																																							
Waist (Inch)	24-28	28-32	32-36	36-40	40-44																																																																																							
<b>Unisex T</b>	Bias decision to chest size - If between sizes, downsizing is preferred																																																																																											
	Small	Medium	Large	X-Large	XX-Large																																																																																							
Chest (Inch)	34-38	38-42	42-46	46-50	50-56																																																																																							
Waist (Inch)	30-34	34-38	38-42	42-46	46-52																																																																																							
<b>Women's Shorts</b>	Bias decision to hip size - If between sizes, downsizing is preferred																																																																																											
	Small	Medium	Large	X-Large	XX-Large																																																																																							
Waist (Inch)	24-28	28-32	32-36	36-40	40-44																																																																																							
Hip (Inch)	34 - 38 <sup>1/2</sup>	38 - 42 <sup>1/2</sup>	42 <sup>1/2</sup> - 46 <sup>1/2</sup>	46 <sup>1/2</sup> - 50 <sup>1/2</sup>	50 <sup>1/2</sup> - 54 <sup>1/2</sup>																																																																																							
<b>Men's Shorts</b>	Bias decision to hip size - If between sizes, downsizing is preferred																																																																																											
	Small	Medium	Large	X-Large	XX-Large																																																																																							
Waist (Inch)	30-34	34-38	38-42	42-46	46-52																																																																																							
Hip (Inch)	35 - 39 <sup>1/2</sup>	39 <sup>1/2</sup> - 43 <sup>1/2</sup>	43 <sup>1/2</sup> - 47 <sup>1/2</sup>	47 <sup>1/2</sup> - 51 <sup>1/2</sup>	51 <sup>1/2</sup> - 55 <sup>1/2</sup>																																																																																							
<b>Wound Care products required</b> <small>(Please Check Mark Size Requested)</small>				<table border="1"> <thead> <tr> <th>Drainage Required</th> <th>Max Units per Wound per Month</th> <th>Frequency of Change Daily unless Spec.</th> <th colspan="8">Wound Number (✓ for supplies)</th> </tr> <tr> <th></th> <th></th> <th></th> <th>Wound #1</th> <th>Wound #2</th> <th>Wound #3</th> <th>Wound #4</th> <th>Wound #5</th> <th>Wound #6</th> <th>Wound #7</th> <th>Wound #8</th> </tr> </thead> <tbody> <tr> <td>Hidrawear Garment (Circle style and size on sizing chart above)</td> <td>-</td> <td>**</td> <td>✓</td> <td>✓</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Hidrawear Dressing Pad Super Absorbent (5x3) <small>(recommended with garment)</small></td> <td>Mod/Heavy</td> <td>30</td> <td>✓</td> <td>✓</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Carboflex: (4x4) (3x6) (6x8)</td> <td>Mod/Heavy</td> <td>30</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Xeroform: (2x2) (4x4) (5x9)</td> <td>-</td> <td>30</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Telfa Non-Adherent Pad: (2x2) (4x4) (4x5) (4x8) Other _____</td> <td>-</td> <td>30</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Gelling Fiber Super Absorbent: (2x2) (4x4) (6x6) (8x8)</td> <td>Mod/Heavy</td> <td>30</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>	Drainage Required	Max Units per Wound per Month	Frequency of Change Daily unless Spec.	Wound Number (✓ for supplies)											Wound #1	Wound #2	Wound #3	Wound #4	Wound #5	Wound #6	Wound #7	Wound #8	Hidrawear Garment (Circle style and size on sizing chart above)	-	**	✓	✓							Hidrawear Dressing Pad Super Absorbent (5x3) <small>(recommended with garment)</small>	Mod/Heavy	30	✓	✓							Carboflex: (4x4) (3x6) (6x8)	Mod/Heavy	30									Xeroform: (2x2) (4x4) (5x9)	-	30									Telfa Non-Adherent Pad: (2x2) (4x4) (4x5) (4x8) Other _____	-	30									Gelling Fiber Super Absorbent: (2x2) (4x4) (6x6) (8x8)	Mod/Heavy	30								
Drainage Required	Max Units per Wound per Month	Frequency of Change Daily unless Spec.	Wound Number (✓ for supplies)																																																																																									
			Wound #1	Wound #2	Wound #3	Wound #4	Wound #5	Wound #6	Wound #7	Wound #8																																																																																		
Hidrawear Garment (Circle style and size on sizing chart above)	-	**	✓	✓																																																																																								
Hidrawear Dressing Pad Super Absorbent (5x3) <small>(recommended with garment)</small>	Mod/Heavy	30	✓	✓																																																																																								
Carboflex: (4x4) (3x6) (6x8)	Mod/Heavy	30																																																																																										
Xeroform: (2x2) (4x4) (5x9)	-	30																																																																																										
Telfa Non-Adherent Pad: (2x2) (4x4) (4x5) (4x8) Other _____	-	30																																																																																										
Gelling Fiber Super Absorbent: (2x2) (4x4) (6x6) (8x8)	Mod/Heavy	30																																																																																										
<b>Practitioner Signature:</b>		<b>Signed Date:</b> 11-2-2022																																																																																										

### Size Matters:

For best user experience please ensure you take your measurements prior to selecting a size. If in between sizes, downsizing is typically recommended. If you are unsure of your size, please take the time to contact our customer care team - call toll-free (888) 711 2014 or email hidrawear@rotech.com

For further support please contact us on  
**Call: (888) 711 2014**  
**Email: hidrawear@rotech.com**

**Rotech WoundCare Complete** is an authorised distributor for HidraWear

<b>Clinic Name</b>	
Clinic Phone	Clinic Fax

<b>Patient Name</b>
---------------------

<b>Order Date:</b>
--------------------





Have patient's wound/s ever been debrided?	YES	NO
--	-----	----

<b>FAX # (866) 220-8963</b>	
Fax or email* with Patient Information Sheet <a href="mailto:orders@rotech.com">orders@rotech.com</a> Phone # (888) 711-2014	
AE Name	AE ID #

Is patient currently seen by Home Health?	YES	NO
---	-----	----

ICD-10	L73.2-Hidradentis suppurativa	L98.499-Non-pressure chronic ulcer of skin	L98.419-Non-pressure chronic ulcer of skin / Buttocks
--------	-------------------------------	--	---

WOUND INFORMATION	ICD-10 (See ICD-10 legend above)	Wound Thickness (HydraWear dressing requires Full)		LOCATION	Wound Size (cm) (L x W x D)	Drainage (Exudate) HydraWear Dressing requires M - H			
		Partial	Full			N	L	M	H
Wound 1:		Partial	Full			N	L	M	H
Wound 2:		Partial	Full			N	L	M	H
Wound 3:		Partial	Full			N	L	M	H
Wound 4:		Partial	Full			N	L	M	H
Wound 5:		Partial	Full			N	L	M	H
Wound 6:		Partial	Full			N	L	M	H
Wound 7:		Partial	Full			N	L	M	H
Wound 8:		Partial	Full			N	L	M	H

Garment Sizing	 <p><b>Crop Top</b></p>	Bias decision to bust size - If between sizes, downsizing is preferred						
		Small	Medium	Large	X-Large	XX-Large		
		Bust (Inch)	32-36	36-40	40-44	44-48	48-52	
	Waist (Inch)	24-28	28-32	32-36	36-40	40-44		
	 <p><b>Unisex T</b></p>	Bias decision to chest size - If between sizes, downsizing is preferred						
		Small	Medium	Large	X-Large	XX-Large		
		Chest (Inch)	34-38	38-42	42-46	46-50	50-56	
	Waist (Inch)	30-34	34-38	38-42	42-46	46-52		
	 <p><b>Women's Shorts</b></p>	Bias decision to hip size - If between sizes, downsizing is preferred						
		Small	Medium	Large	X-Large	XX-Large		
		Waist (Inch)	24-28	28-32	32-36	36-40	40-44	
		Hip (Inch)	34 - 38 <sup>1/2</sup>	38 <sup>1/2</sup> - 42 <sup>1/2</sup>	42 <sup>1/2</sup> - 46 <sup>1/2</sup>	46 <sup>1/2</sup> - 50 <sup>1/2</sup>	50 <sup>1/2</sup> - 54 <sup>1/2</sup>	
	 <p><b>Men's Shorts</b></p>	Bias decision to hip size - If between sizes, downsizing is preferred						
		Small	Medium	Large	X-Large	XX-Large		
		Waist (Inch)	30-34	34-38	38-42	42-46	46-52	
	Hip (Inch)	35 - 39 <sup>1/2</sup>	39 <sup>1/2</sup> - 43 <sup>1/2</sup>	43 <sup>1/2</sup> - 47 <sup>1/2</sup>	47 <sup>1/2</sup> - 51 <sup>1/2</sup>	51 <sup>1/2</sup> - 55 <sup>1/2</sup>		

\*\* One Garment provided with first order

HydraWear is the first and only Hidradentis Suppurativa (HS) specific wound dressing system, intended for use by people with wounds in difficult to dress areas of the body such as the armpit, buttocks & groin. It is comprised of a reusable washable garment along with ergonomically designed single-use wound pads that come with a patented dressing fastener. HydraWear is clinically proven to assure secure dressing retention without the use of adhesives

DRESSING SUPPLIES	Wound Care products required <i>(Please Check Mark Size Requested)</i>	Drainage Required	Max Units per Wound per Month	Frequency of Change Daily unless Spec.	Wound Number (✓ for supplies)													
					Wound #1		Wound #2		Wound #3		Wound #4							
					#	#	#	#	#	#	#	#						
	Hidrawear Garment (Circle style and size on sizing chart above)	-	**															
	Hidrawear Dressing Pad Super Absorbent (5x3) (recommended with garment)	Mod/Heavy	30															
	Carboflex: (4x4) (3x6) (6x8)	Mod/Heavy	30															
	Xeroform: (2x2) (4x4) (5x9)	-	30															
	Telfa Non-Adherent Pad: (2x2) (4x4) (4x5) (4x8) Other_____	-	30															
	Gelling Fiber Super Absorbent: (2x2) (4x4) (6x6) (8x8)	Mod/Heavy	30															

✓	Practitioner Name	Practitioner Name	Practitioner Name	Practitioner Name	Practitioner Name	Practitioner Name	<b>Treating Practitioner MUST check mark their respective box</b>
	NPI#	NPI#	NPI#	NPI#	NPI#	NPI#	

SIGNATURE	* I request that payment of my insurance benefits for any supplies be made to Rotech Healthcare Inc. I am responsible for any balance due that is not covered by my insurance. I understand any product received in my home cannot be returned if opened. I authorize any holder of my medical information to release to Rotech Healthcare Inc. any information needed to determine benefits payable for these supplies or services. I acknowledge receipt and understanding of the patient rights that Rotech Healthcare Inc. provides to all patients.					
	Patient Signature*:					Date:
	I attest by my signature that it is my intention for the prescription to remain valid until the diagnosis described is resolved or otherwise directed by the signer. The requested supplies are medically necessary and the wound(s) has/have been debrided and/or surgically created or modified. I have instructed the patient on how to use the supplies being requested.					Duration of Treatment will be 90 Days unless specified:
<b>X Practitioner Signature:</b>					Signed Date:	

\* If emailing, please use encrypted/secure email. Communications via unencrypted email are not secure and there is a possibility PHI may be misdirected or intercepted.

# Frequently Asked Questions

## What is needed to process an order?

The following information is a general checklist of the items required on a standard written order (Rotech form or any written order) by most insurance plans, including Medicare.

### 1 Qualified Wound Is Present

Medicare covers surgical/wound dressings when a qualifying wound is present. CMS defines a qualified wounds as either of the following:

- A wound caused by, or treated by a surgical procedure
- After debridement of the wound, regardless of the debridement technique

Document debridement TYPE used to remove devitalized or necrotic tissue from wound (examples given are not all-inclusive):

**Surgical:** sharp instruments or laser

**Mechanical:** irrigation or wet-to-dry

**Chemical:** topical application of enzymatic agents

**Autolytic:** application of an occlusive dressing to open wound

### 2 Wound Information

- Type of wound(s) or ICD-10 Codes
- Location of each wound
- Size of wound(s) in cm (L x W x D)
- Amount of exudate

### 3 Treatment Plan

- Type of dressing to use
- Size of the dressings
- Amount to be used at one time
- Frequency of the dressing change
- Expected duration (up to 90 days)

### 4 Provider's Information

- Provider's name and NPI
- If unable to obtain signature when order is placed, please note where signature request should be sent (if it is different than clinic)

### 5 Patient's Demographic Info / Face Sheet

Having this information right away allows us to begin verifying insurance benefits and helps avoid any delays in shipment.

## How easy is it to order and have supplies shipped to patients?

- Rotech's one-page order form takes just 45 seconds to complete. (...yes we timed it!)
- As an added bonus, Rotech is an option when ordering via iHeal, Wound Expert & Intellicure.

## How quickly are supplies delivered?

- With the exception of a PO Box address, orders are shipped and delivered within 24-48 hours

## What if Rotech is not contracted with a patient's insurance plan?

- In most cases, Insurance is verified within an hour of receiving the order, and once approved, supplies ship the same day.
- If Rotech is not contracted with a patient's plan, we will do our very best to find a company that is, notifying you along the way!

## Will Rotech substitute products when I request brand specific products?

- Brand specific products will not be substituted for a lesser product.
- If there is a reimbursement issue or a product is not available, our team will confirm any changes with you prior to shipping supplies.

## Our staff does not have time to deal with unnecessary phone calls. Is this something we should expect?

- Not only does our team value your time but we respect it. You are unlikely to receive any calls from us unless we are returning a call to you.
- Our team will fax order updates to confirm shipment and make you aware of any delays.