



Redefining Patient led Wound Care



Hidradenitis Suppurativa Specific Wound Dressing System

HidraWear is the worlds first **HS specific** wound care product designed to make managing HS drainage quick, easy, and completely adhesive free.

It is a unique **Wearable Wound Care** system for the management of wounds in difficult to dress areas, such as the under-arm region, or buttocks, groin & upper thigh. This **NEW** dressing retention system represents a brand new category in wound care. It combines a reusable garment with a super absorbent dressing, ensuring **adhesive free dressing retention** through our patented Stick-tite technology.


The Garments

Each garment utilizes our patented **HidraHex** technology to enable patients to easily change their dressing day or night.

HidraWear HidraSorb Dressing

8 x 13 cm Super Absorbent dressing locks away moisture from the wound while protecting clothing and bedlinen from exudate strike through.


Stick-tite Hook & Loop retention technology secures the dressing in place without use of any adhesives.




Women's Crop Top




Unisex T-Shirt



Men's Boxer Briefs



Women's Briefs



Retention Dressing System includes a washable garment and non adhesive super absorbent dressings

HidraWear Garments Incorporate:



HidraHex

The hexagon pattern is designed to minimize leaks and strike through while remaining fully breathable & promoting wicking.



Easy On

HidraWear Garments are specifically designed to be easy to get on and off , with the women's underarm crop top incorporating a front closing clasp.



SeamSense

HidraWear garments have outward facing seams. is ensures soft fabric is against your delicate skin and reduces irritation. On the shorts and briefs the seams are on the reverse and provide smooth and frictionless movement.



AeroSilk

All HidraWear garments are made of a premium micro modal and elastane blend fabric that is soft and breathable.



TrueFit

Adjustable straps in the women's underarm crop top enable a tailored fit.

PRIVATE PATIENTS & DPS CARD HOLDERS

Patient led Wound Care

If your patient is on the Drug Payment Scheme, HydraWear products can be prescribed and dispensed through community pharmacy.

Recommended Provision

A starter pack is recommended for the initial prescription with dressing refills as required thereafter.

1. Fill out prescription form
2. Patient brings script to pharmacy
3. Pharmacy orders from Uniphar
(See table for GMS codes)
4. Pharmacy receives order & dispenses to patient

ORANMORE PRIMARY CARE CENTRE
MAIN ST
ORANMORE,
Co. GALWAY H91 V207

TEL: 091 798 1112
FAX: 091 798 1212
EMAIL: prescriptions@OPCC.ie

PRIVATE PRESCRIPTION

PLEASE FILL OUT ALL INFORMATION REQUIRED

PATIENTS NAME:
Ms. Katherine Healy

ADDRESS:
123 Hyde Park, Oranmore, Co Galway

HidraWear Starter Pack, Small, Crop Top
1 Box (60) HidraWear dressings per month for 6 months

SIGNATURE OF PRESCRIBER: Deirdre Kenny

DATE: 01 Apr 2022

DOCTORS NAME: Dr Dierdre Kenny

ADDRESS: Oranmore Primary Care Centre, Oranmore Co Galway

REGISTRATION NUMBER: 42884



Women's Crop Top

Women's Briefs

Men's Boxer Briefs

Unisex T-Shirt

ORDERING INFORMATION FOR PHARMACIST

Covered by PCRS	Crop Top Starter Pack*	Product Code	GMS Code	Ordering Code
	Size Small	HIDRAWEAR SM - STARTERPKAXG-01	86433	787293
	Size Medium	HIDRAWEAR MED - STARTERPKAXG-02	86433	787294
	Size Large	HIDRAWEAR LGE - STARTERPKAXG-03	86433	787295
	Size X-Large	HIDRAWEAR X/LGE - STARTERPKAXG-04	86433	787296
	Size XX-Large	HIDRAWEAR XX/LGE - STARTERPKAXG-05	86433	787297
HidraWear Dressings		Product Code	GMS Code	Ordering Code
13cm x 8cm dressings. (60/box)		BHD001	86432	786565

Starter Pack Contains

- ✓ 2x Dressing retention garments (wear one, wash one)
- ✓ 1x Box of 60 dressings (approx. 1 month supply for 2x exuding wounds)

MEDICAL CARD PATIENTS

Patient led Wound Care

The pink stock order form should be used by all participating GP's to obtain supplies of HydraWear for their medical card patients

Pink Stock Order Forms must contain the following:

- ✓ Pharmacy name/address
- ✓ Practice name/address (on original copy)
- ✓ Signature on the form
- ✓ Computer sequence number in the space provided
- ✓ Date on which the stock order form is issued
- ✓ Clearly indicate in columns 1 and 2 the size and quantity of item(s) required

Pink Stock Order Forms can be ordered from your local community healthcare organisation.

GENERAL MEDICAL SERVICES

TO: PHARMACY ADDRESS ABC Pharmacy
O'Connell St.
Dranmore
Co. Galway

DOCTOR'S NAME ADDRESS SIGNATURE Dr. D. Kenny
Dranmore Primary Care
Dranmore, Co. Galway
Deirdre Kenny

PLEASE SUPPLY ME WITH THE UNITS OF STOCK LISTED IN COLS. 1 AND 2 BELOW

MONTH 04 YEAR 22

STOCK REQUIRED	Quantity Required	Quantity Received	MOH SUPPLIES Drug Code	Computer Units Supplied	This column is for official use
1 <u>HydraWear Starter Pack, Small, Crop Top</u>	<u>1</u>				
2 <u>HydraWear Dressing, 60 Box</u>	<u>3</u>				
3					
4					
5					
6					
7					

RECEIVED THE STOCK ITEMS LISTED IN COLS. 1 AND 2 AND THE QUANTITIES SHOWN IN COL. 3

DOCTOR'S SIGNATURE Deirdre Kenny

DATE 04 / 04 / 2022

FOR OFFICIAL USE BY LOCAL HEALTH OFFICE

321123



1. Fill out pink stock order form & send to pharmacy
2. Pharmacy orders from Uniphar (See table for GMS codes)
3. Pharmacy receives order & sends to GP practice
4. Patient returns to clinic to collect

The declaration at the foot of each order form regarding receipt of stock items **should not be signed and dated until doctors have checked the stock received.**

REQUEST FOR PRESCRIPTION

When to Use

Should you recommend the prescription of HydraWear for the management of your patients wound, please provide the following information to their prescribing physician.

Please arrange for the following items to be:

☐ Ordered and supplied to patient using a pink stock order form (**Medical Card**)

☐ Prescribed (**DPS**)

Item List:

HydraWear Starter Packs contains **2 reusable medical garments** and one pack of **60 wound dressings**

☐ HydraWear Crop Top Starter Pack SIZE **SMALL***
GMS code 86433, Uniphar Ordering Code for pharmacy 787293

☐ HydraWear Crop Top Starter Pack SIZE **X LARGE***
GMS code 86433, Uniphar Ordering Code for pharmacy 787296

☐ HydraWear Crop Top Starter Pack SIZE **MEDIUM***
GMS code 86433, Uniphar Ordering Code for pharmacy 787294

☐ HydraWear Crop Top Starter Pack SIZE **XX LARGE***
GMS code 86433, Uniphar Ordering Code for pharmacy 787297

☐ HydraWear Crop Top Starter Pack SIZE **LARGE***
GMS code 86433, Uniphar Ordering Code for pharmacy 787295

☐ HydraWear Dressings (13cm x 8cm) 60 Pack
GMS Code 86432, Uniphar ordering code for pharmacy 786565

Supply dressings once monthly for ☐ 1 month ☐ 3 months ☐ 6 months

*Starter packs for the **T-Shirt & Briefs retention garments** are not currently covered on GMS. HydraMed will inform when they are available through the primary reimbursement scheme.

Women's Crop Top - Inches

Measurements	Small		Medium		Large		X-Large		XX-Large	
Bust (in)	32	36	36	40	40	44	44	48	48	52
Waist (in)	24	28	28	32	32	36	36	40	40	44

Women's Briefs - Inches

Measurements	Small		Medium		Large		X-Large		XX-Large	
Hip (Widest Part) (in)	34	38.5	38.5	42.5	42.5	46.5	46.5	50.5	50.5	54.5
Waist (in)	24	28	28	32	32	36	36	40	40	44

Men's Boxer Briefs - Inches

Measurements	Small		Medium		Large		X-Large		XX-Large	
Hip (Widest Part) (in)	35	39.5	39.5	43.5	43.5	47.5	47.5	51.5	51.5	55.5
Waist (in)	30	34	34	38	38	42	42	46	46	52

Unisex T-Shirt - Inches

Measurements	Small		Medium		Large		X-Large		XX-Large	
Chest (in)	34	38	38	42	42	46	46	50	50	56
Waist (in)	30	34	34	38	38	42	42	46	46	52

For further support, please contact us at
01-961-7891 or **info@hidrawear.com**

